

# Daily record

---

## Day 1 (July 1): Preparation of Plasmids

### Goal:

- Clone the synthetic INP–PBP and PeIB–PBP sequences into the pET23b plasmid for further transformation into *E. coli* DH5α and BL21 strains.

### Steps:

1. **Plasmid Preparation:** Obtain the pET23b plasmid, linearize it using NdeI and XhoI restriction enzymes.
2. **Gene Insert:** Insert the codon–optimized INP–PBP and PeIB–PBP sequences into the linearized plasmid via ligation.
3. **Transformation into DH5α:**
  - a. Transform the ligated plasmids into *E. coli* DH5α competent cells via heat shock.
  - b. Plate the transformed cells on LB agar plates with Ampicillin (50 µg/mL) for selection.

### Expected Results:

- Colonies should appear after overnight incubation, indicating successful transformation of DH5α with the recombinant plasmids.
- 

## Day 2 (July 2): Transformation and Colony Selection

### Goal:

- Select and verify transformed *E. coli* DH5α colonies containing the recombinant INP–PBP and PeIB–PBP plasmids.

### Steps:

1. **Colony PCR:** Pick individual colonies from the LB agar plates and perform colony PCR to verify the presence of the recombinant plasmids.
2. **PCR Analysis:**
  - a. Use specific primers for INP–PBP and PeIB–PBP to amplify target sequences.
  - b. Run the PCR products on an agarose gel and visualize bands to confirm successful cloning.

### Expected Results:

- PCR results should show bands of expected size for INP–PBP and PeIB–PBP, confirming the presence of the target genes in the selected colonies.
- 

## Day 3 (July 3): Plasmid Extraction and Sequencing

**Goal:**

- Extract plasmids from the positive colonies and send them for sequencing to verify the correct sequence of INP–PBP and PelB–PBP constructs.

**Steps:**

1. **Plasmid Extraction:**
  - a. Use a plasmid extraction kit (e.g., Tiangen) to extract the plasmids from confirmed DH5α colonies.
  - b. Quantify the plasmid concentration using a Nanodrop spectrophotometer.
2. **Sequencing:** Send the extracted plasmids for Sanger sequencing to confirm the correct gene insertion.

**Expected Results:**

- Plasmid yields should be sufficient (>100 ng/μL) for sequencing. Sequencing should confirm the correct integration of INP–PBP and PelB–PBP constructs in the plasmids.
- 

**Day 4 (July 4): Transformation into BL21****Goal:**

- Transform the verified plasmids into *E. coli*/BL21 for protein expression.

**Steps:**

1. **Heat Shock Transformation:**
  - a. Transform the confirmed INP–PBP and PelB–PBP plasmids into BL21 competent cells.
  - b. Incubate the transformed cells on LB agar plates with Ampicillin for selection.
2. **Incubation:** Allow colonies to grow overnight at 37°C.

**Expected Results:**

- Colonies should grow on the Ampicillin plates, indicating successful transformation of BL21 with the recombinant plasmids.
- 

**Day 5 (July 5): Initial Culturing of BL21 Strains****Goal:**

- Begin culturing *E. coli* BL21 strains containing INP–PBP and PelB–PBP to prepare for phosphate adsorption experiments.

**Steps:**

1. **Inoculation:**
  - a. Pick individual colonies from the transformation plates and inoculate them into LB liquid medium containing Ampicillin (50 μg/mL).
2. **Incubation:** Grow the cultures overnight at 37°C, shaking at 250 rpm to ensure optimal growth.

### Expected Results:

- The bacterial cultures should turn cloudy, indicating successful growth of the BL21 strains expressing INP–PBP and PelB–PBP.
- 

## Day 6 (July 6): Plasmid Extraction from BL21 Strains

### Goal:

- Extract plasmids from the transformed BL21 strains to verify the stability of the recombinant plasmids.

### Steps:

1. Plasmid Extraction:
  - a. Use a plasmid extraction kit to isolate the plasmids from the BL21 cultures.
  - b. Measure the plasmid concentration using a Nanodrop spectrophotometer.
2. Plasmid Check: Perform restriction enzyme digestion (NdeI and XhoI) and run on an agarose gel to check the integrity of the plasmids.

### Expected Results:

- Sufficient plasmid yield and correct size bands on the gel should confirm the integrity of the plasmids in BL21 strains.
- 

## Day 7 (July 7): Initial Phosphate Adsorption Test for INP–PBP

### Goal:

- Perform the first phosphate adsorption test for the INP–PBP engineering strain.

### Steps:

1. Culture Preparation:
  - a. Adjust the bacterial culture to an OD600 of 1 in Tris–HCl buffer (pH 7.4).
  - b. Add 10 mg/L KH<sub>2</sub>PO<sub>4</sub> to the solution.
2. Incubation: Incubate the culture at room temperature, shaking at 250 rpm for 3 hours.
3. Phosphate Detection: After incubation, measure the remaining phosphate concentration in the supernatant using the Malachite Green Phosphate Detection Kit.

### Expected Results:

- Phosphate adsorption should be observed, with a decrease in phosphate concentration in the solution.
- 

## Day 8 (July 8): Phosphate Adsorption Test for PelB–PBP

### Goal:

- Perform the same phosphate adsorption experiment as the previous day, but for the PelB–PBP strain.

### Steps:

1. Culture Preparation:
  - a. Adjust the PeIB–PBP bacterial culture to OD600=1.
  - b. Add 10 mg/L KH<sub>2</sub>PO<sub>4</sub> to the solution.
2. Incubation: Incubate at room temperature (250 rpm, 3 hours).
3. Phosphate Detection: Measure the phosphate concentration using the Malachite Green Phosphate Detection Kit.

### Expected Results:

- Phosphate adsorption should be evident, allowing for comparison with the results of INP–PBP.
- 

## Day 9 (July 9): Initial Data Analysis

### Goal:

- Analyze the phosphate adsorption data from both INP–PBP and PeIB–PBP strains.

### Steps:

1. Data Entry:
  - a. Enter the phosphate concentration data into GraphPad Prism for analysis.
2. Comparison: Compare the adsorption capacities of both strains based on the change in phosphate concentration.
3. Statistical Test: Perform a t–test to check if there are significant differences between the adsorption capacities of the two strains.

### Expected Results:

- A preliminary understanding of which strain shows better adsorption efficiency under standard conditions.
- 

## Day 10 (July 10): Optimization of Adsorption (pH Testing for INP–PBP)

### Goal:

- Begin optimizing the phosphate adsorption for INP–PBP by testing its adsorption efficiency at different pH levels.

### Steps:

1. Buffer Preparation:
  - a. Prepare buffers of pH 3, 5, 7, 8, and 10 (using acetic acid for pH 3 and 5, PBS for pH 7, and Tris–HCl for pH 8 and 10).
2. Adsorption Experiment:
  - a. Adjust the OD600 of the bacterial culture to 1 and suspend it in the different buffers.
  - b. Add 10 mg/L KH<sub>2</sub>PO<sub>4</sub> and incubate at room temperature, shaking at 250 rpm for 3 hours.

3. **Phosphate Detection:** Measure phosphate concentrations after incubation using the Malachite Green Kit.

**Expected Results:**

- Identify the pH level at which INP–PBP exhibits the highest phosphate adsorption.
- 

**Day 11 (July 11): Optimization of Adsorption (pH Testing for PeIB–PBP)**

**Goal:**

- Perform the same pH optimization experiment for PeIB–PBP, similar to the INP–PBP strain.

**Steps:**

1. **Buffer Preparation:**
  - a. Prepare buffer solutions at different pH levels (3, 5, 7, 8, 10) using appropriate buffers (acetic acid for pH 3 and 5, PBS for pH 7, Tris–HCl for pH 8 and 10).
2. **Adsorption Experiment:**
  - a. Adjust the OD600 of PeIB–PBP culture to 1.
  - b. Add 10 mg/L KH<sub>2</sub>PO<sub>4</sub> to the buffer solutions and incubate for 3 hours at 250 rpm at room temperature.
3. **Phosphate Detection:** Measure the phosphate concentrations in the supernatant after incubation using the Malachite Green Phosphate Detection Kit.

**Expected Results:**

- Compare the phosphate adsorption capacities at different pH levels to identify the optimal pH for PeIB–PBP.
- 

**Day 12 (July 12): Optimization of Adsorption (Temperature Testing for INP–PBP)**

**Goal:**

- Optimize the adsorption efficiency of the INP–PBP strain by testing at different temperatures.

**Steps:**

1. **Temperature Adjustment:**
  - a. Prepare Tris–HCl buffer (pH 7.4) and adjust the reaction temperature to 25°C, 35°C, and 45°C using temperature–controlled shakers.
2. **Adsorption Experiment:**
  - a. Adjust the OD600 of the bacterial culture to 1, suspend it in the buffer, and add 10 mg/L KH<sub>2</sub>PO<sub>4</sub>.
  - b. Incubate the culture for 3 hours at each temperature (shaking at 250 rpm).
3. **Phosphate Detection:** Measure phosphate concentration in the supernatant using the Malachite Green Phosphate Detection Kit.

**Expected Results:**

- Determine the optimal temperature for phosphate adsorption in INP–PBP. Higher temperatures may enhance adsorption but could also cause denaturation.
- 

## Day 13 (July 13): Optimization of Adsorption (Temperature Testing for PelB–PBP)

### Goal:

- Perform the same temperature optimization experiment for PelB–PBP, as done with INP–PBP.

### Steps:

1. **Temperature Setup:**
  - a. Prepare Tris–HCl buffer (pH 7.4) and adjust the temperature to 25°C, 35°C, and 45°C.
2. **Adsorption Experiment:**
  - a. Adjust the OD600 of PelB–PBP culture to 1, add 10 mg/L KH<sub>2</sub>PO<sub>4</sub>, and incubate for 3 hours at each temperature (250 rpm).
3. **Phosphate Detection:** Use the Malachite Green Kit to measure the remaining phosphate in the solution.

### Expected Results:

- Identify the temperature that provides the highest phosphate adsorption for PelB–PBP.
- 

## Day 14 (July 14): Initial Phosphate Desorption Study for INP–PBP

### Goal:

- Evaluate the phosphate desorption capacity of the INP–PBP strain under different pH and temperature conditions.

### Steps:

1. **Culture Preparation:**
  - a. Culture INP–PBP strains overnight and adjust to OD600=1 in Tris–HCl buffer (pH 7.5).
  - b. Saturate the bacteria with phosphate by incubating in 150 mg/L KH<sub>2</sub>PO<sub>4</sub> for 3 hours.
2. **Desorption Experiment:**
  - a. Wash the bacteria twice with buffer (10 mM Tris–HCl, 1 mM MgCl<sub>2</sub>, pH 7.5) to remove unbound phosphate.
  - b. Incubate the phosphate–saturated bacteria in desorption buffers of varying pH (3, 5, 7, 8, 10) and temperatures (25°C, 35°C, 45°C) for 3 hours at 250 rpm.
3. **Phosphate Detection:** Use the Malachite Green Kit to measure the phosphate concentration in the supernatant after desorption.

### Expected Results:

- Observe how pH and temperature affect phosphate release from the INP–PBP bacteria. Higher temperatures and more extreme pH values (e.g., pH 3) may promote desorption.
- 

## Day 15 (July 15): Phosphate Desorption Study for PelB–PBP

### Goal:

- Perform the same phosphate desorption study for the PelB–PBP strain, similar to INP–PBP.

### Steps:

1. **Culture Preparation:**
  - a. Adjust the PelB–PBP culture to OD<sub>600</sub>=1 and saturate with 150 mg/L KH<sub>2</sub>PO<sub>4</sub> for 3 hours.
  - b. Wash the bacteria to remove unbound phosphate.
2. **Desorption Experiment:**
  - a. Incubate the phosphate–saturated cells in desorption buffers at different pH levels (3, 5, 7, 8, 10) and temperatures (25°C, 35°C, 45°C).
3. **Phosphate Detection:** Measure phosphate release with the Malachite Green Kit.

### Expected Results:

- Evaluate the effect of different desorption conditions on PelB–PBP. Compare results with those from the INP–PBP strain.
- 

## Day 16 (July 16): Isolation of Cytoplasmic and Membrane Components of INP–PBP

### Goal:

- Separate the cytoplasmic and membrane components of INP–PBP engineering bacteria and analyze their respective phosphate binding capacities.

### Steps:

1. **Bacterial Culture:** Grow 50 mL of INP–PBP bacteria overnight and harvest the cells via centrifugation at 10,000 rpm for 10 minutes.
2. **Cell Lysis:**
  - a. Re–suspend the pellet in 10 mL PBS and lyse the cells using ultrasonication (150W, 1 second on, 3 seconds off, for 20 minutes).
3. **Component Separation:**
  - a. Centrifuge the lysate at 5,000 rpm for 10 minutes to remove debris, and then at 39,000 rpm for 1 hour to separate the membrane fraction (pellet) from the cytoplasmic fraction (supernatant).
  - b. Collect both fractions separately for further testing.
4. **Phosphate Adsorption Test:** Test the adsorption capacity of the cytoplasmic and membrane components using Malachite Green Kit.

### Expected Results:

- Determine whether the **membrane component** or **cytoplasmic component** contributes more to phosphate adsorption.
- 

## Day 17 (July 17): Isolation of Soluble and Periplasmic Components of PelB–PBP

### Goal:

- Perform similar component isolation for PelB–PBP, focusing on **soluble** and **periplasmic** fractions.

### Steps:

1. **Bacterial Culture:** Grow 50 mL of PelB–PBP bacteria overnight and harvest cells by centrifugation.
2. **Periplasmic Extraction:**
  - a. Re-suspend the bacterial pellet in **TES buffer** and incubate for 10 minutes.
  - b. Centrifuge at 10,000×g and collect the **periplasmic component** (supernatant).
3. **Soluble Component Extraction:**
  - a. Lyse the bacteria using **sonication** to collect the **soluble cytoplasmic component**.
4. **Phosphate Adsorption Test:** Test the phosphate adsorption ability of each component.

### Expected Results:

- Compare phosphate adsorption between the **periplasmic** and **soluble components** of PelB–PBP. Expect the periplasmic fraction to have higher adsorption, as the PelB signal directs proteins to this region.
- 

## Day 18 (July 18): Data Analysis for Adsorption and Desorption Tests

### Goal:

- Analyze data from the **adsorption** and **desorption** experiments for both INP–PBP and PelB–PBP strains.

### Steps:

1. **Data Entry:** Enter all phosphate concentration data from adsorption and desorption experiments into **GraphPad Prism**.

### Statistical Tests:

1.
  - a. Use **one-way ANOVA** to determine the significance of pH and temperature effects on adsorption/desorption.
  - b. Perform **post hoc Tukey's tests** to compare specific groups.

**Graph Plotting:** Plot graphs for **phosphate adsorption** as a

1. function of **pH** and **temperature** for both strains.

### Expected Results:

- Significant differences between different conditions (pH and temperature), with optimal conditions clearly identified.
- 

## Day 19 (July 19): Replication of Key Experiments

### Goal:

- Repeat key adsorption and desorption experiments to verify reproducibility of results under optimal conditions.

### Steps:

1. **Adsorption:** Perform additional adsorption experiments for both strains at pH 7 and 35°C (optimal conditions identified).
2. **Desorption:** Repeat desorption at pH 7.5 and 45°C for both strains.
3. **Phosphate Detection:** Measure phosphate adsorption/desorption using Malachite Green Kit.

### Expected Results:

- Reproducibility of the previously identified optimal conditions.
- 

## Day 20 (July 20): Final Analysis and Experimental Adjustments

### Goal:

- Conduct a final analysis of all experiments and make any necessary adjustments for ongoing or future experiments.

### Steps:

1. **Final Data Analysis:** Complete the statistical analysis and confirm the key findings.
2. **Troubleshooting:** If certain experiments (e.g., desorption) did not yield consistent results, plan for **troubleshooting steps** or method adjustments.
3. **Planning:** Prepare for the next set of experiments based on the data gathered.

### Expected Results:

- Finalized understanding of the optimal conditions for phosphate adsorption and desorption, with clear next steps for any further experiments.
- 

## Day 21 (July 21): Optimization of Phosphate Adsorption and Desorption for INP-PBP

### Goal:

- Fine-tune and verify the phosphate adsorption and desorption properties of INP-PBP under optimal conditions.

### Steps:

1. **Culture Preparation:**
  - a. Grow INP-PBP strains in LB + Amp (50 µg/mL) at 37°C, 250 rpm overnight.

- b. Adjust the culture to OD<sub>600</sub>=1 using Tris-HCl buffer (pH 7.5).
2. Adsorption Experiment:
  - a. Add 10 mg/L KH<sub>2</sub>PO<sub>4</sub> to the culture and incubate for 3 hours at the optimized temperature (35°C).
3. Desorption Experiment:
  - a. After adsorption, perform the desorption step by suspending phosphate-saturated bacteria in buffer at 45°C and pH 7.5 for 3 hours.
4. Phosphate Measurement: Use Malachite Green Phosphate Detection Kit to measure phosphate levels after adsorption and desorption.

#### Expected Results:

- Final verification of phosphate adsorption and desorption capacities at optimized conditions, confirming pH 7.5 and 35–45°C as ideal for adsorption and desorption.
- 

### Day 22 (July 22): Optimization of Phosphate Adsorption and Desorption for PelB-PBP

#### Goal:

- Repeat and verify the phosphate adsorption and desorption experiments for PelB-PBP at optimal conditions.

#### Steps:

1. Culture Preparation: Grow PelB-PBP strains overnight at 37°C, 250 rpm in LB + Amp (50 µg/mL).
2. Adsorption: Adjust the culture to OD<sub>600</sub>=1 and suspend in Tris-HCl buffer (pH 7.5). Add 10 mg/L KH<sub>2</sub>PO<sub>4</sub> and incubate at 35°C for 3 hours.
3. Desorption: Transfer the phosphate-saturated bacteria to a desorption buffer (pH 7.5) and incubate at 45°C for 3 hours.
4. Phosphate Measurement: Measure phosphate adsorption and desorption using the Malachite Green Kit.

#### Expected Results:

- Similar to INP-PBP, optimal adsorption and desorption conditions should be confirmed at pH 7.5 and temperatures between 35°C – 45°C.
- 

### Day 23 (July 23): Immobilization of INP-PBP on Sepharose Beads

#### Goal:

- Perform immobilization of INP-PBP onto NHS-Activated Sepharose beads to investigate the possibility of creating a reusable phosphate-capturing system.

#### Steps:

1. Immobilization Setup:
  - a. Prepare fresh NHS-Activated Sepharose 4 Fast Flow beads by washing them with 1 mM HCl.

- b. Mix INP–PBP enzyme solution (20 mL) with the washed NHS beads and incubate for 16 hours at 4°C, 30 rpm.
2. **Bead Washing:** Wash the immobilized beads three times with Tris–HCl (pH 7.5) to remove unbound protein.
3. **Phosphate Adsorption Test:**
  - a. Test the adsorption capacity of the INP–PBP immobilized beads by incubating them in a 10 mg/L KH<sub>2</sub>PO<sub>4</sub> solution at room temperature (25°C) for 10 minutes.
4. **Phosphate Measurement:** Use Malachite Green Kit to measure the phosphate concentration in the supernatant.

#### **Expected Results:**

- The immobilized INP–PBP should retain its phosphate adsorption capacity, enabling the beads to act as a reusable phosphate–capturing system.
- 

### **Day 24 (July 24): Immobilization of PelB–PBP on Sepharose Beads**

#### **Goal:**

- Perform the same immobilization experiment for PelB–PBP on Sepharose beads as conducted for INP–PBP.

#### **Steps:**

1. **Immobilization:** Mix PelB–PBP enzyme solution (20 mL) with washed NHS–Activated Sepharose beads and incubate for 16 hours at 4°C, 30 rpm.
2. **Washing:** Wash the immobilized beads three times with Tris–HCl (pH 7.5).
3. **Phosphate Adsorption Test:**
  - a. Incubate PelB–PBP beads in 10 mg/L KH<sub>2</sub>PO<sub>4</sub> at 25°C for 10 minutes to test phosphate adsorption.
4. **Phosphate Measurement:** Measure phosphate concentration in the supernatant using the Malachite Green Kit.

#### **Expected Results:**

- PelB–PBP immobilized beads should demonstrate effective phosphate adsorption, providing a comparison with the INP–PBP system.
- 

### **Day 25 (July 25): Phosphate Desorption Test for Immobilized INP–PBP**

#### **Goal:**

- Test the ability of the INP–PBP immobilized beads to release phosphate under different desorption conditions.

#### **Steps:**

1. **Saturation and Washing:**
  - a. Saturate INP–PBP beads with phosphate by incubating them in 150 mg/L KH<sub>2</sub>PO<sub>4</sub> at 25°C.

- b. Wash the beads three times with Tris-HCl buffer (pH 7.5) to remove unbound phosphate.
2. Desorption Experiment:
  - a. Incubate the beads in buffers of varying pH (3, 5, 7, 8, 10) and temperatures (25°C, 35°C, 45°C) for 10 minutes to assess phosphate release.
3. Phosphate Measurement: Measure the amount of phosphate released from the beads using the Malachite Green Kit.

**Expected Results:**

- Optimal desorption conditions should be identified, with higher temperatures and extreme pH values likely promoting phosphate release.
- 

## Day 26 (July 26): Phosphate Desorption Test for Immobilized PeIB-PBP

**Goal:**

- Perform the same phosphate desorption test for PeIB-PBP immobilized beads as conducted for INP-PBP.

**Steps:**

1. Saturation and Washing:
  - a. Saturate PeIB-PBP beads with phosphate and wash them to remove unbound phosphate.
2. Desorption: Incubate the beads in different pH and temperature conditions to test phosphate release.
3. Phosphate Measurement: Use the Malachite Green Kit to quantify phosphate released from the beads.

**Expected Results:**

- PeIB-PBP beads should show effective phosphate release, allowing for comparison with INP-PBP beads.
- 

## Day 27 (July 27): Comparison of Adsorption Efficiency Between Free and Immobilized PBPs

**Goal:**

- Compare the phosphate adsorption efficiency of free INP-PBP and PeIB-PBP in solution with the immobilized PBPs on Sepharose beads.

**Steps:**

1. Adsorption Comparison:
  - a. Test the adsorption of free INP-PBP and PeIB-PBP in solution using 10 mg/L KH<sub>2</sub>PO<sub>4</sub> under optimal conditions (pH 7.5, 35°C).
  - b. Repeat the adsorption test with immobilized PBPs under the same conditions.
2. Phosphate Measurement: Measure the remaining phosphate in solution after adsorption for both free and immobilized forms using the Malachite Green Kit.

### Expected Results:

- Immobilized PBPs may exhibit slightly reduced adsorption compared to free PBPs due to steric hindrance, but immobilization offers the advantage of reusability.
- 

## Day 28 (July 28): Statistical Analysis and Data Visualization

### Goal:

- Conduct comprehensive statistical analysis and visualize the results from adsorption, desorption, and immobilization experiments.

### Steps:

1. **Data Entry:** Organize all phosphate concentration data from adsorption, desorption, and immobilization tests.
2. **Statistical Tests:**
  - a. Perform ANOVA to assess the impact of different conditions (pH, temperature) on adsorption and desorption efficiency.
  - b. Use Tukey's post hoc test for pairwise comparison of groups.
3. **Data Visualization:** Plot adsorption/desorption curves and bar graphs comparing free and immobilized PBP efficiency.

### Expected Results:

- Significant differences between conditions should be evident, with clear identification of optimal conditions for both free and immobilized PBP systems.
- 

## \*\*Day

29 (July 29)\*\*: *\*Experimental Review and Troubleshooting\**

### Goal:

- Review the experimental results, troubleshoot any issues, and plan for potential reruns or method adjustments.

### Steps:

1. **Review Results:** Analyze the overall findings from the phosphate adsorption, desorption, and immobilization experiments.
2. **Troubleshooting:** If any data appears inconsistent or unclear, plan for **replicating** experiments under controlled conditions.
3. **Method Adjustments:** Identify any procedural weaknesses and make adjustments for future experiments if necessary.

### Expected Results:

- Clarification of any experimental inconsistencies and a clear plan for potential reruns.
- 

## Day 30 (July 30): Final Report Preparation and Conclusion

**Goal:**

- Summarize the experimental findings, prepare a report, and finalize the results for presentation to the evaluation committee.

**Steps:****1. Summary of Findings:**

- a. Prepare a comprehensive summary of the phosphate adsorption and desorption capacities of both INP–PBP and PeIB–PBP, along with the results from the immobilization experiments.

**2. Conclusions:**

- a. Conclude which conditions (pH, temperature) are optimal for phosphate capture and release for each strain.
- b. Discuss the advantages and limitations of **immobilized PBP** systems.

**3. Report Preparation:** Compile all data, analysis, and conclusions into a formal report for presentation.**Expected Results:**

- A clear, well–organized report that presents the experimental data and conclusions, highlighting the success of the INP–PBP and PeIB–PBP engineering systems for phosphate adsorption.